

OR10X1 Antibody



PACO02473

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR10X1

Uniprot

Q8NGY0/Q8NGX7

Synonyms:

Olfactory receptor 10X1

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 10X1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COLO cells using Olfactory receptor 10X1 Polyclonal Antibody.

