

APOOL Antibody



PACO02481

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

APOOL

Uniprot

Q6UXV4

Synonyms:

APOOL; CXorf33; FAM121A; Apolipoprotein O-like; Protein FAM121A


Immunogen:

Synthesized peptide derived from the C-terminal region of human ApoOL.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of SH-SY5Y 293T 3T3 lysis using ApoOL antibody.