

# OR2AG2/OR2AG1 Antibody



PACO02487

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:****Gene ID:**

OR2AG2/OR2AG1

**Uniprot**

A6NM03/Q9H205

**Synonyms:**

OR2AG2; OR2AG2P; Olfactory receptor 2AG2; OR2AG1; OR2AG3; Olfactory receptor 2AG1; HT3; Olfactory receptor 2AG3; Olfactory receptor OR11-79

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 2AG1/2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 293 cells using Olfactory receptor 2AG1/2 Polyclonal Antibody.

