

# CPNE8 Antibody



PACO02488

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

CPNE8

**Uniprot**

Q86YQ8

**Synonyms:**

CPNE8; Copine-8; Copine VIII

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human Copine 8.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 3T3 cells using Copine 8 Polyclonal Antibody.

