

CSNK1D Antibody



PACO02493

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

CSNK1D

Uniprot

P48730

Synonyms:

CSNK1D; HCKID; Casein kinase I isoform delta; CKI-delta; CKId; Tau-protein kinase CSNK1D

Immunogen:

Synthesized peptide derived from the Internal region of human Casein Kinase Idelta.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of brain cells using Casein Kinase Idelta Polyclonal Antibody.

