

# PDIK1L Antibody



PACO02516

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

PDIK1L

**Uniprot**

Q8N165

**Synonyms:**

PDIK1L; CLIK1L; Serine/threonine-protein kinase PDIK1L; PDLIM1-interacting kinase 1-like

**Immunogen:**

Synthesized peptide derived from the Internal region of human PDIK1L.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of K562 cells using PDIK1L Polyclonal Antibody.

