

ALS2CR11 Antibody



PACO02525

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

C2CD6

Uniprot

Q53TS8

Synonyms:

ALS2CR11; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 11 protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human ALS2CR11.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 3T3 cells using ALS2CR11 Polyclonal Antibody.

