

BTLA Antibody



PACO02529

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

BTLA

Uniprot

Q7Z6A9

Synonyms:

BTLA; B- and T-lymphocyte attenuator; B- and T-lymphocyte-associated protein; CD antigen CD272

Immunogen:

Synthesized peptide derived from the Internal region of human CD272.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of K562 cells using CD272 Polyclonal Antibody.

