

# GUSBP1 Antibody



PACO02531

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

GUSBP1

**Uniprot**

Q15486

**Synonyms:**

GUSBP1; SMA3; Putative beta-glucuronidase-like protein SMA3

**Immunogen:**

Synthesized peptide derived from the Internal region of human GUSBP1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of K562 cells using GUSBP1 Polyclonal Antibody.

