

# Phospho-ADRBK1 (S685) Antibody



PACO02538

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

GRK2

**Uniprot**

P25098

**Synonyms:**

ADRBK1; BARK; BARK1; GRK2; Beta-adrenergic receptor kinase 1; Beta-ARK-1; G-protein coupled receptor kinase 2

**Immunogen:**

Synthesized peptide derived from human GRK 2 around the phosphorylation site of S685.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of JK cells using Phospho-GRK 2 (S685) Polyclonal Antibody.

