

FAM84B Antibody



PACO02545

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

FAM84B

Uniprot

Q96KN1

Synonyms:

FAM84B; BCMP101; NSE2; Protein FAM84B; Breast cancer membrane protein 101; Protein NSE2

Immunogen:

Synthesized peptide derived from the N-terminal region of human NSE2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of Jurkat cells using NSE2 Polyclonal Antibody.

