

TMEM30B Antibody



PACO02561

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TMEM30B

Uniprot

Q3MIR4

Synonyms:

TMEM30B; CDC50B; Cell cycle control protein 50B; Transmembrane protein 30B

Immunogen:

Synthesized peptide derived from the Internal region of human Cdc50B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of K562 cells using Cdc50B Polyclonal Antibody.

