

GPR113 Antibody



PACO02572

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ADGRF3

Uniprot

Q8IZF5

Synonyms:

GPR113; PGR23; Probable G-protein coupled receptor 113; G-protein coupled receptor PGR23

Immunogen:

Synthesized peptide derived from the N-terminal region of human GPR113.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HeLa cells using GPR113 Polyclonal Antibody.

