

# KCNJ16 Antibody



PACO02573

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

KCNJ16

**Uniprot**

Q9NPI9

**Synonyms:**

KCNJ16; Inward rectifier potassium channel 16; Inward rectifier K(+) channel Kir5.1; Potassium channel; inwardly rectifying subfamily J member 16

**Immunogen:**

Synthesized peptide derived from mouse KIR5.1 around the non-phosphorylation site of S416.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HeLa cells using KIR5.1 Polyclonal Antibody.

