

KCNV2 Antibody



PACO02578

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

KCNV2

Uniprot

Q8TDN2

Synonyms:

KCNV2; Potassium voltage-gated channel subfamily V member 2; Voltage-gated potassium channel subunit Kv8.2

Immunogen:

Synthesized peptide derived from the Internal region of human KV8.2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of MCF-7 cells using KV8.2 Polyclonal Antibody.

