

PWWP2B Antibody



PACO02581

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

PWWP2B

Uniprot

Q6NUJ5

Synonyms:

PWWP2B; PWWP2; PWWP domain-containing protein 2B

Immunogen:

Synthesized peptide derived from the Internal region of human PWWP2B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of LOVO cells using PWWP2B Polyclonal Antibody.

