

S100Z Antibody



PACO02582

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

S100Z

Uniprot

Q8WXG8

Synonyms:

S100Z; Protein S100-Z; S100 calcium-binding protein Z

Immunogen:

Synthesized peptide derived from the N-terminal region of human S-100Z.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COLO cells using S-100Z Polyclonal Antibody.

