

# OR5M3 Antibody



PACO02694

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR5M3

**Uniprot**

Q8NGP4

**Synonyms:**

OR5M3; Olfactory receptor 5M3; Olfactory receptor OR11-191

**Immunogen:**


Synthesized peptide derived from the Internal region of human Olfactory receptor 5M3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HEPG2 cells using Olfactory receptor 5M3 Polyclonal Antibody.