

OR8G5 Antibody



PACO02696

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR8G5

Uniprot

Q8NG78

Synonyms:

OR8G5; OR8G5P; OR8G6; Olfactory receptor 8G5; Olfactory receptor 8G6; Olfactory receptor OR11-298


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 8G5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Olfactory receptor 8G5 Polyclonal Antibody.