

OR9Q1 Antibody



PACO02698

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR9Q1

Uniprot

Q8NGQ5

Synonyms:

OR9Q1; Olfactory receptor 9Q1


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 9Q1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Olfactory receptor 9Q1 Polyclonal Antibody.