

GPR114 Antibody



PACO02705

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

ADGRG5

Uniprot

Q8IZF4

Synonyms:

GPR114; PGR27; Probable G-protein coupled receptor 114; G-protein coupled receptor PGR27


Immunogen:

Synthesized peptide derived from the Internal region of human GPR114.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using GPR114 Polyclonal Antibody.