

# NACAD Antibody



PACO02754

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

NACAD

**Uniprot**

O15069

**Synonyms:**

NACAD; KIAA0363; NAC-alpha domain-containing protein 1

**Immunogen:**


Synthesized peptide derived from the Internal region of human NACAD.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat cells using NACAD Polyclonal Antibody.