

GPS2 Antibody



PACO02789

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

GPS2

Uniprot

Q13227

Synonyms:

GPS2; G protein pathway suppressor 2; GPS-2


Immunogen:

Synthesized peptide derived from the N-terminal region of human GPS2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAW264.7 cells using GPS2 Polyclonal Antibody.