

Phospho-GYS1 (S645) Antibody



PACO02805

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

GYS1

Uniprot

P13807

Synonyms:

GYS1; GYS; Glycogen [starch] synthase; muscle


Immunogen:

Synthesized peptide derived from human Glycogen synthase 1 around the phosphorylation site of S645.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Phospho-Glycogen Synthase 1 (S645) Polyclonal Antibody.