

# ICAM2 Antibody



PACO02829

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

ICAM2

**Uniprot**

P13598

**Synonyms:**

ICAM2; Intercellular adhesion molecule 2; ICAM-2; CD antigen CD102

**Immunogen:**

Synthesized peptide derived from the Internal region of human ICAM-2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 3T3 cells using ICAM-2 Polyclonal Antibody.

