

# OR2W3 Antibody



PACO02835

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

OR2W3

**Uniprot**

Q7Z3T1

**Synonyms:**

OR2W3; OR2W3P; OR2W8P; Olfactory receptor 2W3; Olfactory receptor 2W8; Olfactory receptor OR1-49

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 2W3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 3T3 cells using Olfactory receptor 2W3 Polyclonal Antibody.

