

# OR5B12 Antibody



PACO02862

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR5B12

**Uniprot**

Q96R08

**Synonyms:**

OR5B12; OR5B12P; OR5B16; Olfactory receptor 5B12; Olfactory receptor 5B16; Olfactory receptor OR11-241

**Immunogen:**


Synthesized peptide derived from the Internal region of human Olfactory receptor 5B12.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of A549 cells using Olfactory receptor 5B12 Polyclonal Antibody.