

PACO02864

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

BSX

Uniprot

Q3C1V8

Synonyms:

BSX; BSX1; Brain-specific homeobox protein homolog

Immunogen:

Synthesized peptide derived from the C-terminal region of human BSX.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293T cells using BSX Polyclonal Antibody.

