

OR5H15 Antibody



PACO02883

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR5H15

Uniprot

A6NDH6

Synonyms:

OR5H15; Olfactory receptor 5H15


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 5H15.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO cells using Olfactory receptor 5H15 Polyclonal Antibody.