

# MAGEA5 Antibody



PACO02887

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

MAGEA5

**Uniprot**

P43359

**Synonyms:**

MAGEA5; MAGE5; Melanoma-associated antigen 5; Cancer/testis antigen 1.5; CT1.5; MAGE-5 antigen

**Immunogen:**

Synthesized peptide derived from the Internal region of human MAGE-A5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 3T3 cells using MAGE-A5 Polyclonal Antibody.

