

# SLC39A4 Antibody



PACO02995

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

SLC39A4

**Uniprot**

Q6P5W5

**Synonyms:**

SLC39A4; ZIP4; Zinc transporter ZIP4; Solute carrier family 39 member 4; Zrt- and Irt-like protein 4; ZIP-4

**Immunogen:**

Synthesized peptide derived from the Internal region of human ZIP4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of Colo cells using ZIP4 Polyclonal Antibody.

