

ENAH Antibody



PACO02997

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

ENAH

Uniprot

Q8N8S7

Synonyms:

ENAH; MENA; Protein enabled homolog

Immunogen:

Synthesized peptide derived from the C-terminal region of human Mena.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using Mena Polyclonal Antibody.

