

Phospho-PRKCD (T507) Antibody



PACO03004

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

PRKCD

Uniprot

Q05655

Synonyms:

PRKCD; Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta


Immunogen:

Synthesized peptide derived from human PKC delta around the phosphorylation site of T507.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 3T3 cells using Phospho-PKC delta (T507) Polyclonal Antibody.