

MDM1 Antibody



PACO03029

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MDM1

Uniprot

Q8TC05

Synonyms:

MDM1; Nuclear protein MDM1

Immunogen:

Synthesized peptide derived from the C-terminal region of human MDM1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 3T3 cells using MDM1 Polyclonal Antibody.

