

ZNF695 Antibody



PACO03031

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ZNF695

Uniprot

Q8IW36

Synonyms:

ZNF695; Zinc finger protein 695; Zinc finger protein SBZF3

Immunogen:

Synthesized peptide derived from the Internal region of human ZNF695.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293T cells using ZNF695 Polyclonal Antibody.

