

RGS12 Antibody



PACO03065

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

RGS12

Uniprot

O14924

Synonyms:

RGS12; Regulator of G-protein signaling 12; RGS12

Immunogen:

Synthesized peptide derived from the N-terminal region of human RGS12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A