

# Phospho-TOP2A (T1343) Antibody



PACO03114

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

TOP2A

**Uniprot**

P11388

**Synonyms:**

TOP2A; TOP2; DNA topoisomerase 2-alpha; DNA topoisomerase II; alpha isozyme

**Immunogen:**

Synthesized peptide derived from human Topo IIalpha around the phosphorylation site of T1343.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A