

MMTAG2 Antibody



PACO03135

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MMTAG2

Uniprot

Q9BU76

Synonyms:

MMTAG2; C1orf35; Multiple myeloma tumor-associated protein 2; hMMTAG2

Immunogen:

Synthesized peptide derived from the N-terminal region of human MMTAG2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A