

SLC52A2 Antibody



PACO03137

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

SLC52A2

Uniprot

Q9HAB3

Synonyms:

SLC52A2; GPR172A; PAR1; RFT3; Solute carrier family 52; riboflavin transporter, member 2; Porcine endogenous retrovirus A receptor 1; PERV-A receptor 1; Protein GPR172A; Riboflavin transporter 3; hRFT3

Immunogen:

Synthesized peptide derived from the Internal region of human GPR172A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A