

Phospho-ATRIP (S68) Antibody



PACO03151

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ATRIP

Uniprot

Q8WXE1

Synonyms:

ATRIP; AGS1; ATR-interacting protein; ATM and Rad3-related-interacting protein

Immunogen:

Synthesized peptide derived from human ATRIP around the phosphorylation site of S68.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A