

Phospho-STAM2 (Y192) Antibody



PACO03170

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

STAM2

Uniprot

O75886

Synonyms:

STAM2; HBP; Signal transducing adapter molecule 2; STAM-2; Hrs-binding protein

Immunogen:

Synthesized peptide derived from human STAM2 around the phosphorylation site of Y192.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A