

# Phospho-IRS1 (Y896) Antibody



PACO03243

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

IRS1

**Uniprot**

P35568

**Synonyms:**

IRS1; Insulin receptor substrate 1; IRS-1

**Immunogen:**

Synthesized peptide derived from human IRS-1 around the phosphorylation site of Y896.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A