

Phospho-IRS1 (S323) Antibody



PACO03244

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

IRS1

Uniprot

P35568

Synonyms:

IRS1; Insulin receptor substrate 1; IRS-1

Immunogen:

Synthesized peptide derived from human IRS-1 around the phosphorylation site of S323.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A