

# PRKAA2 Antibody



PACO03305

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PRKAA2

**Uniprot**

P54646

**Synonyms:**

5'-AMP-activated protein kinase catalytic subunit alpha-2

**Immunogen:**

Recombinant Protein.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---

N/A

N/A