

PRKAB2 Antibody



PACO03306

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

Protein Background:

Gene ID:

PRKAB2

Uniprot

O43741

Synonyms:

5'-AMP-activated protein kinase subunit beta-2

Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A