

# OPRK1 Antibody



PACO03311

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

OPRK1

**Uniprot**

P41145

**Synonyms:**

OPRK1; OPRK; Kappa-type opioid receptor; K-OR-1; KOR-1

**Immunogen:**

Synthesized peptide derived from mouse KOR-1 around the non-phosphorylation site of S369.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A