

CST8 Antibody



PACO03317

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

CST8

Uniprot

O60676

Synonyms:

CST8; CRES; Cystatin-8; Cystatin-related epididymal spermatogenic protein

Immunogen:

Synthesized peptide derived from the Internal region of human Cystatin 8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A