

CDH10 Antibody



PACO03318

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:**Gene ID:**

CDH10

Uniprot

Q9Y6N8

Synonyms:

CDH10; Cadherin-10; T2-cadherin

Immunogen:

Synthesized peptide derived from the Internal region of human K-cadherin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A