

# EIF1B Antibody



PACO03332

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

EIF1B

**Uniprot**

O60739

**Synonyms:**

EIF1B; Eukaryotic translation initiation factor 1b; eIF1b; Protein translation factor SUI1 homolog GC20

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human eIF1B.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A