

SNAPC5 Antibody



PACO03336

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

SNAPC5

Uniprot

O75971

Synonyms:

SNAPC5; SNAP19; snRNA-activating protein complex subunit 5; SNAPc subunit 5; Small nuclear RNA-activating complex polypeptide 5; snRNA-activating protein complex 19 kDa subunit; SNAPc 19 kDa subunit

Immunogen:

Synthesized peptide derived from the N-terminal region of human SNAPC 19.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A